

5-50L Rotary Evaporator

Item Number: KE-10



Introduction

Efficiently separate low-boiling solvents with the 5-50L Rotary Evaporator. Ideal for chemical labs, it offers precise and safe evaporating processes.

[Learn More](#)

Model	KR-1005	KR1010	KR1020	KR1050
Evaporating Flask	5L	10L	20L	50L
Flask Size	ø50mm	ø50/125mm	ø125mm	ø125mm
Evaporating Speed	2L/H by H2O	3.2L/H by H2O	5L/H by H2O	9L/H by H2O
Recycling Flask	3L	5L	10L	20L
Condenser	Vertical	Diagonal		
Vacuum Degree	≤2kPa/h			
Vacuum Sealing	Teflon+Viton Sealing			
Speed Control	Knop Stepless Speed Regulation			
Motor Power	250W			
Rotating Speed	20-140rpm	20-130rpm		20-110rpm
Lifting Method	Electric			
Lifting Distance	150mm	160mm	190mm	180mm