

Desktop Fast Autoclave Sterilizer 16L / 24L

Item Number: KT-B02



Introduction

The desktop fast steam sterilizer is a compact and reliable device used for rapid sterilization of medical, pharmaceutical, and research items.

[Learn More](#)

Model	KT-T16J	KT-T24J
Sterilization chamber volume	16L ϕ 250×360mm	24L ϕ 270×410mm
Fixed work stress	0.22Mpa	
Rated working temperature	134°C	
Sterilization temperature selection	105-134°C	
Timing time selection range	0-99min	
Heat uniformity of sterilization room	$\leq \pm 1^\circ\text{C}$	
Power/supply voltage	AC220V.50Hz/2.7KW	
Overall dimensions (mm)	545×328×423	560×350×500
Packing size (mm)	605×390×480	620×410×570
Gross weight/net weight	17/15Kg	20/18Kg