

No Demolding Lab Infrared Press Mold

Item Number: PMI



Introduction

Effortlessly test your samples with no demolding required using our lab infrared press mold. Enjoy high transmittance and customizable sizes for your convenience.

[Learn More](#)

Model	PMI-A	PMI-B
Material	Carbide YT15	
Indenter hardness	HRC85-HRC90	
Sample size	φ13 mm	φ7 mm
Cavity depth	10 mm	5 mm
External dimensions	76×50×70 mm	76×30×70 mm
Weight(Kg)	0.75	0.35

Other sizes can be customized