

# Xrf Boric Acid Lab Powder Pellet Pressing Mold

Item Number: PMXB



## Introduction

Get accurate results with our XRF Boric Acid lab Powder Pellet Pressing Mold. Perfect for preparing samples for X-ray fluorescence spectrometry. Custom sizes available.

[Learn More](#)

Model	PMXB-A
Material	Cr12MoV
Indenter hardness	HRC60-HRC62
Sample size	φ32 / φ40 mm
Cavity depth	45 mm
External dimensions	φ73×133 mm
Weight(Kg)	3.2

**Other sizes can be customized**